



碧云天生物技术/Beyotime Biotechnology  
订货热线: 400-1683301或800-8283301  
订货e-mail: order@beyotime.com  
技术咨询: info@beyotime.com  
网址: http://www.beyotime.com

## Histone H3 Mouse Monoclonal Antibody

产品编号	产品名称	包装
AF0009	Histone H3 Mouse Monoclonal Antibody	100μl

### 产品简介:

来源	应用	交叉反应性	分子量
Mouse	WB, ICC, IF	H, M, R	17KDa

WB, Western blot; IP, Immunoprecipitation; IF, Immunofluorescence; IHC, Immunohistochemistry; ICC, Immunocytochemistry; FC, Flow Cytometry; ELISA, Enzyme-linked Immunosorbent Assay; ChIP, Chromatin Immunoprecipitation Assay.

H, Human; M, Mouse; R, Rat; C, Chicken; Cw, Cow; Dg, Dog; Gp, Guinea pig; Hm, Hamster; Hr, Horse; Mk, Monkey; Pg, Pig;

Rb, Rabbit; S, Sheep; Z, Zebrafish; All, all species expected.

➤ 配套提供了Western一抗稀释液，可以用于Western检测或其它适当用途时的一抗稀释。

➤ 建议抗体使用时的稀释比如下(实际使用时需根据抗原水平的高低作适当调整):

WB	IP	IF	IHC	ICC	FC	ELISA	ChIP
1:1000	-	1:100	-	1:50	-	-	-

### ➤ 抗体详细信息如下:

About this Antibody	
Name	Histone H3 Mouse Monoclonal Antibody
Category	Monoclonal antibody (mAb); Primary antibody
Isotype	IgG2b
Purification	Affinity purification
About the Immunogen	
Immunogen	synthetic peptide (KLH-coupled)
Gene ID	3020/126961/3021
SwissProt	P68431/P68433/Q6LEDO
Synonyms	H3 histone antibody, HIST1H3A antibody, Histone cluster 1, H3a antibody
Category	Chromatin Regulation / Nuclear Function
Background	The nucleosome, made up of DNA wound around eight core histone proteins (two each of H2A, H2B, H3, and H4), is the primary building block of chromatin. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling.

### 包装清单:

产品编号	产品名称	包装
AF0009	Histone H3 Mouse Monoclonal Antibody	100μl
AZ100	Western一抗稀释液	100ml
—	说明书	1份

### 保存条件:

Histone H3 Mouse Monoclonal Antibody -20°C保存, Western一抗稀释液-20°C或4°C保存, 一年有效。Western一抗稀释液优先推荐4°C保存, 长期不使用可以考虑-20°C保存, 但冻融可能会导致出现轻微的浑浊和少量不溶物。

### 注意事项:

➤ 如果本抗体用于Western blot (WB)、免疫荧光(IF)、免疫细胞化学(ICC)等实验, 请注意回收使用过的稀释抗体。回收的抗体通

常至少可以重复使用5-10次。稀释后的抗体，包括已经使用过的稀释抗体，请4°C保存。

- 回收后重复使用的抗体，使用方法同新鲜稀释的抗体。如果在重复使用过程中发现抗体出现轻微混浊现象，可以10,000g离心1-3分钟，取上清用于后续检测。如果回收的抗体出现明显的絮状物或长霉长菌等情况，则可以考虑废弃该抗体。
- 提供的Western一抗稀释液也可以用于免疫荧光(IF)、免疫组化(IHC)、免疫细胞化学(ICC)等适当用途。如果希望获得最佳的检测效果，请考虑使用上述检测专用的一抗稀释液。
- 本产品仅限于专业人员的科学的研究用，不得用于临床诊断或治疗，不得用于食品或药品，不得存放于普通住宅内。
- 为了您的安全和健康，请穿实验服并戴一次性手套操作。

## 使用说明：

请根据抗体的实际用途选择相应的使用方法。

### 1. Western检测：

- a. 按照推荐的稀释比例用碧云天提供的Western一抗稀释液稀释抗体。
- b. 把经过封闭的蛋白膜与稀释好的一抗4°C缓慢摇动过夜或室温缓慢摇动2小时，确保稀释的抗体至少能在摇动的瞬间覆盖蛋白膜。
- c. 回收稀释的一抗，4°C保存以备下次继续使用。
- d. 按照Western的实验步骤进行后续的洗涤、二抗孵育、洗涤和检测等操作。具体操作可以参考如下网页：  
<http://www.beyotime.com/support/western.htm>

### 2. 免疫染色：

可以使用碧云天生产的免疫染色一抗稀释液(P0103)稀释抗体，使用后注意回收稀释好的一抗，具体操作可以参考如下网页：  
<http://www.beyotime.com/support/immunol-staining.htm>

### 3. 其它实验操作请自行参考适当的protocol进行。

### 4. 代表性图片：




Fig. 1. Western blot analysis on HeLa (1), PC12 (2) and F9 (3) cell lysates using Histone H3 mouse monoclonal antibody (AF0009).




Fig. 2. Immunofluorescent staining of HeLa cells using Histone H3 mouse monoclonal antibody (AF0009).




Fig. 3. Immunofluorescent staining of F9 cells using Histone H3 mouse monoclonal antibody (AF0009).

## 相关产品：

产品编号	产品名称	包装
P0006	Bradford蛋白浓度测定试剂盒	1000次

P0010	BCA蛋白浓度测定试剂盒(增强型)	500次
P0012	BCA蛋白浓度测定试剂盒	500次
P0012A	SDS-PAGE凝胶配制试剂盒	1盒
P0012AC	SDS-PAGE凝胶快速配制试剂盒	1盒
P0013	Western及IP细胞裂解液	100ml
P0013B	RIPA裂解液(强)	100ml
P0014B	SDS-PAGE电泳液	10×1L
P0015	SDS-PAGE蛋白上样缓冲液(5X)	2ml
P0018	BeyoECL Plus(超敏ECL化学发光试剂盒)	共100ml
P0018A	BeyoECL Star(特超敏ECL化学发光试剂盒)	共100ml
P0020	显影定影试剂盒	各1升
P0021B	Western转膜液	10×1L
P0023A	Western一抗稀释液	100ml
P0023B	Western封闭液	100ml
P0023C6	Western洗涤液(10X)	10×100ml
P0023D	Western二抗稀释液	100ml
P0025	Western一抗二抗去除液	250ml
P0252	QuickBlock™ Western封闭液	100ml
P0256	QuickBlock™ Western一抗稀释液	100ml
P0258	QuickBlock™ Western二抗稀释液	100ml
P0066	预染蛋白质分子量标准	200μl
P0068	彩色预染蛋白质分子量标准	200μl
P0071	BeyoColor™彩色预染蛋白质分子量标准	200μl
P0098	免疫染色固定液	100ml
P0102	免疫染色封闭液	100ml
P0103	免疫染色一抗稀释液	100ml
P0106L	免疫染色洗涤液(10X)	10×100ml
P0108	免疫荧光染色二抗稀释液	100ml
P0110	免疫染色(非荧光)二抗稀释液	100ml
P0260	QuickBlock™免疫染色封闭液	100ml
P0262	QuickBlock™免疫染色一抗稀释液	100ml
P0265	QuickBlock™免疫荧光染色二抗稀释液	100ml
P0267	QuickBlock™免疫组化染色二抗稀释液	100ml
P0126	抗荧光淬灭封片液	5ml
FFP24	PVDF膜(进口分装, 6.6×8.5cm, 0.2μm)	20张/包装
FFP32	PVDF膜(进口原装, 6.6×8.5cm, 0.45μm)	20张/包装
FFP51	转印滤纸(7.5×10cm)	100片/包装
FF057	X-OMAT BT胶片(原柯达, 5×7英寸)	100张/盒
FFC58	压片暗盒(5×7英寸)	1个

## 使用本产品的文献：

1. Dejun Xu, Lingbin Liu, Yongju Zhao, Li Yang, Jianyong Cheng, Rongmao Hua, Zelin Zhang, Qingwang Li . Melatonin protects mouse testes from palmitic acid-induced lipotoxicity by attenuating oxidative stress and DNA damage in a SIRT1-dependent manner J Pineal Res. 2020 Nov;69(4):e12690.
2. Jinliang Guo, Yang Zhang, Huiling Gao, Shengben Li, Zhen-Yu Wang, Chao-Feng Huang . Mutation of HPR1 encoding a component of the THO/TREX complex reduces STOP1 accumulation and aluminium resistance in Arabidopsis thaliana New Phytol. 2020 Oct;228(1):179-193.
3. Nannan Zheng, Xuanqi Liu, Ying Yang, Yifei Liu, Furong Yan, Yiming Zeng, Yunfeng Cheng, Duoqiao Wu, Chengshui Chen, Xiangdong Wang . Regulatory roles of NAT10 in airway epithelial cell function and metabolism in pathological conditions Cell Biol Toxicol. 2023 Aug;39(4):1237-1256.
4. Sun L, Li M, Ma X, Feng H, Song J, Lv C, He Y . Inhibition of HMGB1 reduces rat spinal cord astrocytic swelling and AQP4 expression after oxygen-glucose deprivation and reoxygenation via TLR4 and NF-κB signaling in an IL-6-dependent manner. J NEUROINFLAMM. 2017 Nov 25;14(1):231.
5. Chen KL, Li L, Li CM, Wang YR, Yang FX, Kuang MQ, Wang GL . SIRT7 Regulates Lipopolysaccharide-Induced Inflammatory Injury by Suppressing the NF-κB Signaling Pathway. Oxid Med Cell Longev. 2019 Jun 11 2019:3187972.
6. Zhang Y, Cui W, Zhai Q, Zhang T, Wen X . N-acetylcysteine ameliorates repetitive/stereotypic behavior due to its antioxidant

- properties without activation of the canonical Wnt pathway in a valproic acid-induced rat model of autism. Mol Med Rep. 2017 Aug;16(2):2233-2240.
7. Wang J, Liu X, Yan C, Liu J, Wang S, Hong Y, Gu A, Zhao P . LEF1-AS1, a long-noncoding RNA, promotes malignancy in glioblastoma. ONCOTARGETS THER. 2017 Aug 28;10:4251-4260.
  8. Fang DZ, Wang YP, Liu J, Hui XB, Wang XD, Chen X, Liu D . MicroRNA-129-3p suppresses tumor growth by targeting E2F5 in glioblastoma. EUR REV MED PHARMACO. 2018 Feb;22(4):1044-1050.
  9. Lin X, Liu X, Gong C . Expression of engrailed homeobox 2 regulates the proliferation, migration and invasion of non-small cell lung cancer cells. Oncol Lett. 2018 Jul;16(1):536-542.
  10. Liu X, Xu Q, Mei L, Lei H, Wen Q, Miao J, Huang H, Chen D, Du S, Zhang S, Zhou J, Deng R, Li Y, Li C, Li H . Paeonol attenuates acute lung injury by inhibiting HMGB1 in lipopolysaccharide-induced shock rats. Int Immunopharmacol. 2018 Aug;61:169-177.
  11. Li D, Guo L, Deng B, Li M, Yang T, Yang F, Yang Z . Long non-coding RNA HR1 participates in the expression of SREBP-1c through phosphorylation of the PDK1/AKT/FoxO1 pathway. Mol Med Rep. 2018 Sep;18(3):2850-2856.
  12. Yu Y, Zhang X, Tian H, Zhang Z, Tian Y . Knockdown of long non-coding RNA HOTAIR increases cisplatin sensitivity in ovarian cancer by inhibiting cisplatin-induced autophagy. J BUON. 2018 Sep-Oct;23(5):1396-1401.
  13. Wu X, Xiao J, Zhao C, Zhao C, Han Z, Wang F, Yang Y, Jiang Y, Fang F . Claudin1 promotes the proliferation, invasion and migration of nasopharyngeal carcinoma cells by upregulating the expression and nuclear entry of  $\beta$ -catenin. Exp Ther Med. 2018 Oct;16(4):3445-3451.
  14. Hui Wang, Zhengxia Wang, Qingbin Tang . Reduced Expression of microRNA-199a-3p Is Associated With Vascular Endothelial Cell Injury Induced by Type 2 Diabetes Mellitus Exp Ther Med. 2018 Oct;16(4):3639-3645.
  15. Song Y, Wang R, Li LW, Liu X, Wang YF, Wang QX, Zhang Q . Long non-coding RNA HOTAIR mediates the switching of histone H3 lysine 27 acetylation to methylation to promote epithelial-to-mesenchymal transition in gastric cancer. Int J Oncol. 2019 Jan 54(1):77-86.
  16. Zhang C, Wang T, Wu H, Zhang L, Li K, Wang F, Chen Y, Jin J, Hua D . HEF1 regulates differentiation through the Wnt5a/ $\beta$ -catenin signaling pathway in human gastric cancer. BIOCHEM BIOPH RES CO. 2019 Jan 29 509(1):201-208.
  17. Sun L, Li M, Ma X, Zhang L, Song J, Lv C, He Y . Inhibiting High Mobility Group Box-1 Reduces Early Spinal Cord Edema and Attenuates Astrocyte Activation and Aquaporin-4 Expression after Spinal Cord Injury in Rats. J NEUROTRAUM. 2019 Feb 1;36(3):421-435.
  18. Dang YF, Qiu TX, Song DW, Liu L . PMA-triggered PKC $\epsilon$  activity enhances Nrf2-mediated antiviral response on fish rhabdovirus infection. FISH SHELLFISH IMMUN. 2019 Apr 87:871-878.
  19. Hua Qian, Yan Cao, Junfeng Sun, Jianing Zu, Liang Ma, Haizhou Zhou, Xianling Tang, Yan Li, Haiyang Yu, Mingyu Zhang, Yunlong Bai, Chaoqian Xu, Norito Ishii, Takashi Hashimoto, Xiaoguang Li . Anti-human serum albumin autoantibody may be involved in the pathogenesis of autoimmune bullous skin diseases FASEB J. 2020 Jun;34(6):8574-8595.
  20. Xudong Zhang, Ting Lian, Wenjun Fan, Guangwei Zhang, Ziwei Chen, Xingchun Gou, Rajiv Kumar Jha . Long-Noncoding RNA CASC9 Promotes Progression of Non-Small Cell Lung Cancer by Promoting the Expression of CDC6 Through Binding to HuR Cancer Manag Res. 2020 Sep 25;12:9033-9043.
  21. Qianqian Li, Ningbo Li, Yeting Zeng, Xinrui Wang, Jie Li, Hongying Su, Meiqin Gao, Xiongfei Huang . Nuclear receptor FXR impairs SK-Hep-1 cell migration and invasion by inhibiting the Wnt/ $\beta$ -catenin signaling pathway Oncol Lett. 2020 Nov;20(5):161.
  22. Guangzhao Song, Liyan Li, Ying Yang . MicroRNA-329-3p alleviates high glucose-induced endothelial cell injury via inhibition of the TLR4/TRAFF/NF- $\kappa$ B signaling pathway Exp Ther Med. 2021 Jan;21(1):29.
  23. Danmei Tian, Lingyun Cao, Qingwen Li, Huiyun Huang, Wei Xu, Guodong Chen, Zhaohui Song, Yi He, Xinsheng Yao, Jinshan Tang . Benzannulated 5,5-spiroketal sesquiterpenes from the roots of Angelica Pubescens Bioorg Chem. 2021 Feb;107:104604.
  24. Ling Zhou, Zeming Liu, Sichao Chen, Jing Qiu, Qianqian Li, Shipei Wang, Wei Zhou, Danyang Chen, Guang Yang, Liang Guo . Transcription factor EB-mediated autophagy promotes dermal fibroblast differentiation and collagen production by regulating endoplasmic reticulum stress and autophagy-dependent secretion Int J Mol Med. 2021 Feb;47(2):547-560.
  25. Guan-Yu Wang, Zhi-Lin Luan, Ning-Wei Che, De-Bin Yan, Xiao-Wan Sun, Cong Zhang, Jian Yin . Inhibition of microRNA-129-2-3p protects against refractory temporal lobe epilepsy by regulating GABRA1 Brain Behav. 2021 Jul;11(7):e02195.
  26. Xi Su, Wenqiang Li, Zhen Li, Kang Liu, Meng Song, Minglong Shao, Luxian Lv, Xulu Chang . Chronic developmental exposure to low-dose ([C8mim][PF6]) induces neurotoxicity and behavioural abnormalities in rats Ecotoxicol Environ Saf. 2021 Dec 1;225:112806.
  27. Zhifang Li, Hexiang Xu, Yi Xu, Guanfeng Lu, Qiwei Peng, Jiefang Chen, Rentang Bi, Jianzhuang Li, Shengcai Chen, Hongkai Li, Huijuan Jin, Bo Hu . Morinda officinalis oligosaccharides alleviate depressive-like behaviors in post-stroke rats via suppressing NLRP3 inflammasome to inhibit hippocampal inflammation CNS Neurosci Ther. 2021 Dec;27(12):1570-1586.
  28. Gaohong Liu, Xiaojun Ren, Yousong Li, Han Li . Midkine promotes kidney injury in diabetic kidney disease by increasing neutrophil extracellular traps formation Ann Transl Med. 2022 Jun;10(12):693.
  29. Shiyao Xu, Yan Xiong, Rui Yao, Rong Tian, Zhuqing Meng, Mohamed Y Zaky, Beibei Fu, Dong Guo, Lulu Wang, Feng Lin, Xiaoyuan Lin, Haibo Wu . A Novel ERK2 Degrader Z734 Induces Apoptosis of MCF-7 Cells via the HERC3/p53 Signaling Pathway Molecules. 2022 Jul 6;27(14):4337.
  30. Heng Jiang, Nan Zhang, Huichen Ge, Jiaping Wei, Xin Xu, Xiangkun Meng, Kun Qian, Yang Zheng, Jianjun Wang . S6K1 acts through FOXO to regulate juvenile hormone biosynthesis in the red flour beetle, *Tribolium castaneum* J Insect Physiol. 2022 Jul:140:104405.
  31. Hanna Dai, Xiaoling Chen, Daiwen Chen, Bing Yu, Jun He, Hong Chen, Hui Yan, Ping Zheng, Yuheng Luo, Zhiqing Huang . Effects of dietary l-theanine supplementation on pork quality and muscle fiber type transformation in finishing pigs J Sci Food Agric. 2023 Mar 15;103(4):2106-2115.
  32. Ruiqing Sun, Shuangying Li, Bin Xia, Junxia Zhu . Detection of novel variant and functional study in a Chinese family with nonsyndromic oligodontia Oral Dis. 2023 Jul;29(5):2177-2187.

Version 2024.03.12